

RADAR

A discreet, humane and innovative solution to mouse control.

Designed specifically with high risk customers in mind. RADAR detects, captures, and isolates mice, efficiently, hygienically and humanely.

This innovative system was developed by Rentokil to meet the needs of high risk market segments, where the presence of a single mouse cannot be tolerated. RADAR is therefore ideal for food & pharmaceutical manufacturers, telecommunications, hospitality, educational and health establishments. It provides discreet and effective pest control which is safer than conventional baits, as there is no risk of pesticide contamination.



Key Customer Benefits

- Effective 24/7 detection and eradication of mice within your premises.
- Highly effective protection, allowing use in any sensitive areas, including production environments.
- Extremely discreet & unobtrusive, units are placed at wall floor junctions.
- Innovative, self-contained unit allows service with no disruption to your business.
- Simple, hygienic, easy to clean & more discreet than other pest control alternatives.
- Can be used in conjunction with Rentokil's PestNetOnline monitoring & reporting system.



How does RADAR work?

Our qualified and experienced Rentokil Technicians will identify rodent runs and fix the RADAR units in appropriate locations.


RADAR has two entrances, one at each end of the unit, allowing mice to run through its passageway. When the mouse breaks two consecutive infrared beams, it trips a circuit that immediately closes both entrances. Once the mouse is trapped inside the tightly sealed chamber, the carbon dioxide gas is released. The gas acts extremely rapidly and is widely recognised as a humane method of control.

There is no danger of contamination as the mouse remains completely isolated.

A warning light immediately illuminates to indicate capture.

As the unit is self-contained it can be safely & discreetly removed from your premises by our qualified and experienced Technicians.

On-going contract maintenance will ensure the optimal performance of RADAR, and full compliance with your health & safety obligations.



Technical Specifications

Dimensions (L x D x H)	308mm x 196mm x 60mm
Unit weight (kg)	415g (including full gas canister)
Detection system	Pulsed, reflective dual infrared beam sensors

To arrange a free survey of your premises, or to gain further information on RADAR please contact us.

Call us on
0800 085 7329

Visit us at
rentokil.co.uk

Use biocides safely. Always read the label and product information before use. RADAR contains the active ingredient carbon dioxide.